



LECTEUR DE GLYCEMIE

Glucose meter set Alphacheck Professional

Automatic encoding, result available in just 5 seconds, only a small amount of blood needed (0.5 µl), blood sampling possible at different body locations, Bio-Sensor Technology, 500 memory locations, test strip ejector. Contents: 1 Glucose meter, 10 strips, 10 lancets, 1 lancing device, 2 batteries, etui



Item	Description
1133140	Alphacheck professional

Blood glucose test strips

for Alphacheck professional. PU = 50 pcs / Box



Item	Description
1133141	Test strips

Pen Needles - Droplet®

sterile, compatible with type A insulin pens, PU = 100 pcs / box



Item	Description
1133142	6 mm x 31 G
1133143	8 mm x 31 G
1133144	10 mm x 29 G
1133145	12 mm x 29 G

