



CENTRIFUGEUSE UNIVERSELLE

8 * 20 ml

ECOspin

Table Top Centrifuge is convenient in operation through direct entry of time and speed using water-resistant keypad.

Saving of previous working parameters, even after switching of. It's tough body, made out of ABS thermoplastic makes cleaning easier than ever.

Automatic stop function, when opening the lid, makes ECOSpin save in use.

Rotor : 8 * 20 ml

Power supply AC 230 V / 50,60 Hz

Power Input 135 W

Max. RCF (relative centrifugal force) 2325 g

Max. Speed 4000 min⁻¹

Timer : 0 – 99 min

Min. Speed 100 min⁻¹

Radio Shielding EN 61326

Height x Width x Depth 245 x 240 x 320 mm

Height with open lid 400 mm

Weight without rotor 10 kg

Max. kinetic energy 980 Nm

Noise level <65dB(A)

