



## COLLECTEUR D'INFORMATION DE PERFUSION SK-900 I

**La solution pour la surveillance et le contrôle des appareils à perfusions**

The shape of SK-900I is like a mobile phone, it can be carried by the nurse in the pocket.

We add a special module in the pump, when the pump under working, the module can sent the infusion information to sk-900I wirelessly.

On sk-900I, it can receive the information such as Bed NO., occlusion, infusion rate, low battery, infusion finish, near finish and so on. It can connected with 50 sets working pumps.

Once an alarm is sent form the pump, it will be displayed on the LCD of the data collector terminal.

Very convenient for the nurse to carry, menu operation mode, easy to operate.



**Please contact us for more details**