



DEFIBRILATEUR AED

Défibrillateur automatique externe A-10

L'AED A-10 est un défibrillateur destiné au public.

Simple, instructions parlantes et pictogrammes facilitent son utilisation.

Une analyse rapide et la technologie biphasique assurent une efficacité maximum.

Ne nécessite pas d'entretien

Livré avec sac de transport

Spécifications (voir page suivante)





DEFIBRILATEUR AED



Dimensions 320 (H) x 260 (W) x 115 (D) mm, about 3 kg



Specifications

Defibrillation	Waveform: Biphasic Smartshot (Impedance compensation) Energy: 50-220 J Charging time: max. 15 sec. Operating mode: Semi-Auto.
ECG	Lead: II Impedance: 25 – 175 Ohm Heart Rate: 20 - 300 bpm VF/VT detection Analysis time: 8 – 13 sec. Display speed: 25mm/s
Indication	LCD type Icon type Status LCD Voice prompt Text CPR Indication
Pads	Storage: 28 months Cable length: 1,5 m
Self-Test	Every 24 hours, 1 week, 1 month Power on self-test Test result: Status LCD
Data Backup	SD card IrDA LAN (option)/WLAN (option)/CDMA (option)
Battery	Type: LiMnO ₂ Storage: 5 years Discharge: Min. 200 Monitoring: ca. 10 hours Voltage: 15 V
Regulatory	AAMI DF80, EN 60601 - 1, EN 60601 - 2 2005 AHA/ERC protocol
Dimensions	HxWxD: ca. 320x260x115 mm
Weight	3kg
Environmental Conditions	Operating: 0-43°C Storage condition: -30-65°C

