



DEVELOPEUSE DENTAIRE SCP-100

Open-room operation, small volume, low noise, simple operation, safe, fast, and durable.

Classified as Wall M counted Type and Desk Top Type, and former can save planar space (Standard configuration).

The latter is movable (Additionality optional)

Mute-film developing disposal may be competed when the developer and fixer is at 18-25°C.

One place of film enters per 30 seconds for mull-film developing disposal, and the development may last 2-8 minutes, and this may be integrally adjusted in accordance with the exposure conditions of x-ray machine.

Multiple films of 2×3, 3×4 and 2.7× 5.4cm may be applicable for the treatment of development.

It is convenient to regularity maintain this machine, and no tool is required to replace chemical solution and wash relatable parts. More cost is saved if the machine is used to develop films.



The product is elaborately manufactured for specialized no nasality industry by further optimized design on the basis of similar products produced by advanced technology in the world, thus we make another contribution to the undertaking of no nasality and human kind

Specification:

Nominal-Voltage:	220V±10%
Nominal Frequency:	50Hz
Nominal Input power:	<12W
Film Run Through Time:	2-8 min
Protecting Category:	11
Defend Grade:	IP21
Weight of unit:	7.5kg
Size:	42×19×45(cm)