



GASTROSCOPE - GAMME



CE

	GIS-WP-98	GIS-WP-100	GIS-8	GIS-10
Field of view:	105°	105°	100°	100°
Observation range:	3-100 mm	3-100 mm	3-100 mm	3-100 mm
Viewing direction:	Straight	Straight	Forward	Forward
Distal end diameter:	9.8 mm	8.8 mm	10.5 mm	9.5 mm
Flexible section diameter:	9.0 mm	9.0 mm	11 mm	9.8 mm
Forceps channel diameter:	2.8 mm	2.8 mm	2.8 mm	2.8 mm
Bending capability Up/Down:	200°/90°	200°/90°	180°/90°	200°/90°
Bending capability Right/Left:	100°/100°	100°/100°	120°/210°	100°/100°
Working length:	1,030 mm	1,030 mm	1,030 mm	1,060 mm
Total length:	1,350 mm	1,350 mm	1,350 mm	1,350 mm
Waterproof	Yes	Yes	No	No

