



LAMPE FRONTALE

Specifications

Medical Head Lamp BT - 410

Main body	■ Illumination : 15,000 lx (at 25 cm working distance) ■ Color Temperature : 6,000 Kelvin ■ Max LED operating life : 50,000 hours ■ Weight : 177 g
Battery pack	■ Operation time : 4 hours ■ Recharge : 800 times ■ Output : 3.7V ■ Capacity : 2,200 mAh ■ Type : Li-ion (built-in protected circuit module) ■ Weight : 46 g
Charger	■ Input : AC 100 ~ 240 V , 50 / 60 Hz ■ Output : DC 4.2 V / 750 mA ■ Charging time : 4 hours
Standard	■ EN60601 - 1 , EN60601 - 1 - 2
Headband size	■ 534 ~ 638 mm
Standard Accessory	■ Main body (Power LED) 1 EA ■ Headband 1 EA ■ Battery pack (1 cell) 1 EA ■ Charger 1 EA ■ Operator's manual 1 EA

- high light output of up to 15,000 Lux
- color temperature of 6,000 ° Kelvin
- powerful LED with long life
- Li-Ion battery for 4 hours continuous operation
- Weighs only 223 g with battery

Includes:

Head light, adjustable headband, Li-Ion battery with charger, manual

