



MONITEUR MULTI PARAMETRES J-1 NEW

Tout en un, **ECG, PNI, SaO2, FC et T°** sur un grand affichage **LCD couleur**

- Compact, léger, facilité de manipulation
- Logiciel pour raccordement au PC
- Historique des données et courbes de tendances
- Ecran tactile en couleur
- Mémorisation de 129 heures d'informations
- Coupure automatique de l'alimentation
- Batteries intégrées NiMH avec indicateur de charge

SPECIFICATIONS

ECG

Input:	3-lead ECG
Lead selection:	I; II; III
Gain selection:	x1, x2, x4, x8
Sweep speed:	12.5, 25, 50mm/s
Bandwidth monitor or operation)	0.5~40Hz (for 0.05~100Hz (for diagnosis)
Calibration Signal:	1 mv
Protection:	Against electro-surgical interference and defibrillation
Heart rate	Measure Range: 20~250bpm
	Accuracy: ±1bpm
	Resolution: 1bpm
Alarm mode:	Audible and visual alarm



NIBP

Method:	Oscillometric
Mode:	Manual, Auto, STAT
Measuring Interval in AUTO Mode:	1, 2, 3, 5, 10, 15, 20, 30, 60 (Min)
Measuring Period in STAT Mode:	5 Min
Measuring Type:	SYS, DIA, MEAN
Resolution:	1mmHg
Over-pressure protection:	Yes

Oxygen saturation (SpO2)

Range:	35~100%
Accuracy:	±2% (70~100%), ±3% (40~70%)
Pulse rate	Range: 0-250 bpm
Accuracy:	±1 bpm

Temperature

Measuring range:	0~50 °C
Resolution:	0.1 °C
Accuracy:	± 0.1 °C

Power Option

Ni/MH	rechargeable battery
AC power:	100-240V, 50-60Hz, 6VA

Dimension / Net Weight

151x128x45 mm / 480g