



RESPIRATEUR D'ANESTHESIE 100B

Portable, Multifonctionnel Adultes et enfants O₂ + N₂O Circuit Manuel, ouvert, semi ouvert et fermé

Main Unit:

Mode:Semi-Open, Semi-Closed or Closed System

Gas Supplies:O₂/ N₂O 0.28MPa ~ 0.6Mpa

Flowmeter: O₂/ N₂O two-tube flowmeter, 0 to 10 L/min(with inter-lock protection devices ensuring the Oxygen concentration > 25%)

O₂ Flush:25 to 75L/Min

O₂ deficiency Alarm: Audible alarm for at least 7 seconds

Respiration Loop:Good Airtight,Easy for using, Easy for cleaning and sterilizing

Vaporizer:

Flow Range:0.2 to 10L/Min

Concentration:0 to 5.0%

Enflurane, Isoflurane or Fluothane

