



VENTILATEUR SOINS INTENSIFS 300B

Le 300B est un ventilateur à commande électronique **avec compresseur intégré**

Spécifications

Le respirateur est livré avec humidificateur, chariot roulant et bras articulé.



Respiratory frequency: IPPV: 5~60bpm SIMV:1~12bpm	
Inspiratory flow:	Infant: 4~24L/min(continuous)
	Adult: 10~120L/min(inspiratory phase)
Inspiratory time: 0.2~3s [for IPPV (VC) and SIMV]	
Peak Pressure: 1~9kPa	
Inhaled oxygen concentrateion: 21~100%	
PEEP/CPAP: 0~1.5kPa	
Pressure support: 0~3.5kPa	
Trigger sensitivity: 0.1~1.0kPa	
Ventilation modes:	Adult: IPPV (VC), IPPV(PC), SIMV, SIMV+PSV, CPAP, PSV
	Infant: IPPV (PC), IMV, CPAP continuous flow, pressure limited, time cycle Emergency: IPPV (PC) constant pressure (O2 concentration, inspiratroy pressure key active only)
MV monitoring	0~40L/min
	Upper alarm limit:5~40L/min
	Lower alarm limit:0~20L/min
Airway pressure monitoring: Pmax Pmean, Pplat, Peep	
Ventilation monitoring: Respiratory Frequency, Minute Volume, Tidal Volume	
Alarms: air supply deficiency, oxygen deficiency, hose detachment, high airway pressure, low airway pressure, high MV, low MV, apnea, power failure.	