



# RESPIRATEURS POUR URGENCE



	SIRIO PLUS 3000	SIRIO PLUS	SIRIO S2/T	SAVEQUICK ■	SIRIO BABY 100-200
<b>OPERATION MODES</b>	IPPV (AUT) + AST PCV PSV + APNOEA BACK-UP SIMV CPAP TWIN PAP SIGH NEB	IPPV (AUT) + AST PCV PSV + APNOEA BACK-UP SIMV CPAP MANUAL	IPPV (AUT) + AST PSV CPAP APNOEA BACK-UP	IPPV (AUT) + AST CPAP IPPV (AUT)	CPAP IMV IPPV (AUT) PCV MANUAL
<b>MINUTE VOLUME</b>	FROM 1 TO 40 L/MIN	FROM 1 TO 25 L/MIN	FROM 1 TO 16 L/MIN	FROM 1 TO 16 L/MIN	FIX 8,5 L/MIN
<b>TIDAL VOLUME</b>	FROM 20 TO 3000 ML	FROM 100 TO 1500 ML	FROM 15 TO 3000 ML	FROM 15 TO 3000 ML	
<b>BREATHING RATE</b>	FROM 5 TO 150 BPM	FROM 5 TO 140 BPM	FROM 5 TO 70 BPM	FROM 5 TO 70 BPM	T. INSP: FROM 0,2 TO 3 SEC. T. EXP: FROM 0,4 TO 30 SEC.
<b>INSPIRATORY PAUSE</b>	FROM 0 TO 50% INSPIRATORY TIME	FROM 0 TO 50% INSPIRATORY TIME			
<b>I:E RATIO</b>	FROM 1:4 TO 3:1	FROM 1:4 TO 3:1	1:1,5	1:1,5	
<b>PRESSURE LIMIT</b>	FROM 2 TO 80 CM H <sub>2</sub> O	FROM 2 TO 80 CM H <sub>2</sub> O	FROM 0 TO 50 CM H <sub>2</sub> O	FROM 0 TO 50 CM H <sub>2</sub> O	FROM 0 TO 60 CM H <sub>2</sub> O
<b>PEEP</b>	FROM 0 TO 20 CM H <sub>2</sub> O	FROM 0 TO 20 CM H <sub>2</sub> O	FROM 0 TO 20 CM H <sub>2</sub> O		FROM 0 TO 20 CM H <sub>2</sub> O
<b>TRIGGER</b>	FROM 1 TO 9 MBAR BELOW PEEP LEVEL	FROM 1 TO 9 MBAR BELOW PEEP LEVEL	FROM -1 TO -6 CM H <sub>2</sub> O	FROM -1 TO -6 CM H <sub>2</sub> O	
<b>O<sub>2</sub> MIXER</b>	FROM 21% TO 99%	FROM 21% TO 99%	21%, 60%, 100%	21%, 60%, 100%	FROM 46% TO 100% (BABY 100) FROM 21% TO 100% (BABY 200)
<b>MEASURED PARAMETERS</b>	FLOW MAX-MEAN PEEP INSP/EXP. TIME PAUSE RATE EXP. MIN VOLUME EXP. TIDAL VOLUME	FLOW MAX-MEAN PEEP INSP/EXP. TIME PAUSE RATE EXP. MIN VOLUME EXP. TIDAL VOLUME	MAX PRESS	MAX PRESS	MAX RATE
<b>DISPLAYED WAVEFORMS</b>	PRESSURE VOLUME LOOPS				
<b>DIMENSION AND WEIGHT</b>	mm, 360X330X300 Kgs, 7,8	mm, 360X180X220 Kgs, 4,9	mm, 230X150X130 Kgs, 3,5	mm, 230X150X130 Kgs, 3,5	mm, 280X150X130 Kgs, 4,5

Pour plus d'informations, merci de nous contacter : [info@excotrained.be](mailto:info@excotrained.be)