



TAMBOURS INOX

TAMBOURS POUR LA STERILISATION

Acier inoxydable 18/10,
Couvercle à charnières et fermetures latérales.



Dimensions disponibles (les plus courants)

Dia x Haut mm	Dia x Haut mm	Dia x Haut mm	Dia x Haut mm
125X80	190X190	290X145	340X290
125X100	190X240	290X160	340X340
125X125	190X290	290X190	340X390
150X80	240X120	290X240	390X145
150X100	240X145	290X290	390X165
150X130	240X160	290X340	390X190
150X150	240X190	290X390	390X240
150X190	240X240	340X145	390X290
190X125	240X290	340X160	390X340
190X145	240X340	340X190	390X390
190X165	240X390	340X240	