



# TUBES PRELEVEMENT SANGUIN



## Vacuum blood collection tubes, sterile

PET, cap PP and Bromobutyl, PU = 100 pcs / rack, 12 racks / carton

Item	Description
------	-------------

- |         |  |
|---------|--|
| 1130670 | K3 EDTA for 3 ml blood, violet stopper<br>13 x 75 mm                                   |
| 1130671 | Potassium-Fluoride + NA2 EDTA for 2 ml blood, grey stopper<br>13 x 75 mm               |
| 1130672 | Lithium heparin for 3 ml blood, green stopper<br>13 x 75 mm                            |
| 1130673 | Sodium citrate 3,2 %, for 2,7 ml blood, blue stopper<br>13 x 75 mm                     |
| 1130674 | Coagulation accelerator for 4 ml blood, red stopper<br>13 x 75 mm                      |
| 1130675 | Coagulation accelerator for 6 ml blood, red stopper<br>13 x 100 mm                     |
| 1130676 | Separating gel and coagulation accelerator for 3,5 ml blood, red stopper<br>13 x 75 mm |
| 1130677 | Separating gel and coagulation accelerator for 5 ml blood, red stopper<br>13 x 100 mm  |



## Holder for vacuum tubes

Disposable, PU = 250 pcs / carton

Item	Description
------	-------------

- |         |                 |
|---------|-----------------|
| 1130678 | made of plastic |
|---------|-----------------|



## BD Precisionglide™ - Cannula, green

21G x 1 1/2 (0,8 x 38 mm), PU = 250 pcs / Box, 10 Boxes / carton

Item	Description
------	-------------

- |         |         |
|---------|---------|
| 1130679 | Cannula |
|---------|---------|