



UNITE ORL COMPLETE 920



E.N.T TREATMENT UNIT

4 spray guns + hot air + ventilation system + suction system




CE 0434
ISO 9001
ISO 13485

Designed For Better Life



UNITE ORL COMPLETE 920



STANDARD ACCESSORIES

| | | | |
|---|---------------------------------|----|---|
| 1 | Spray gauge | 10 | Electric light power switch |
| 2 | Suction gauge | 11 | Ventilation system |
| 3 | Hot air | 12 | Suction unit power switch |
| 4 | Spray regulator | 13 | Electric treatment chair power switch |
| 5 | Suction regulator | 14 | Exhaust power switch |
| 6 | Medicine vials shelf | 15 | Main power switch |
| 7 | 4 spray guns | 16 | Electric light with flexible adjustment |
| 8 | Hot air power switch | 17 | Foot control switch socket |
| 9 | Compression system power switch | | |

SPECIFICATIONS

| | | | |
|------------------|--------------------|--------------------|------------------------|
| Dimensions | L81 X W50 X H84 cm | Compressor Pump | 200W |
| Power Supply | AC 110/220V | Max Pressure | 1.1 kg/cm ² |
| Electric Current | 15/10A | Spray Gun | 4 |
| Bulb No. | 1 | Ventilation system | Yes |
| Light Bulb | 4.8V 0.5A | Hot Air | Yes |
| Vacuum Pump | 200W | N.W. | 87 KGS |
| Max Suction | 580 mm/hg | G.W. | 135 KGS |
| Max Volume | 3,000 + 1,000c.c. | Packing Size | L87 X W77 X H135 cm |



UNITE ORL COMPLETE 920

Configuration

Unit with 4 castors,
Main construction: SS304,
Oil-less 200W twin piston **vacuum** pump
Oil-less 200W twin piston **compressor** pump
1pc, 3,000c.c. PC bottle
1pc, 1,000c.c. PC bottle x
1pc, Suction tube
4 pc Spray gun
Ventilation system with Eustachian tube catheter
1set, Mirror warmer
1set, Mini surgeon head light with battery pack 4.8V 0.5A
1pc, Spare liquid medicine bottle
6pcs, Tongue depressor
2 pcs of storage slots
1 pc Instrument tray
1 pc Big round cotton pot
1 pc Small round