



LAMPE D'EXAMEN LED

M 130 LED

Light intensity in Lux at 1 meter distance	50.000
Colour rendering index $R_{a(t)}$ at 4500 Kelvin	95
Focussable light field size (in cm)	17 fixe
Colour temperature (Kelvin)	4500
Electronic light intensity control at the lamphead	50-100%
Temperature increase in head area	0,5 °C
Total power consumption	20 W
Number of LED's	19
Life-span of the LED's	≥ 40.000 h
Working distance (in cm)	70 - 140
Diameter of the light head (in cm)	33
Height adjustment (in cm)	119



Existe en 4 versions

- Montée sur statif roulant à 5 branches
- Montée sur bras articulé fixé sur statif roulant
- Montée sur bras articulé avec fixation murale
- Montée sur bras articulé avec fixation plafond

Option :

- Réglage du focus de 14 à 25 cm