



## AUTOCLAVE ELECTRONIQUE 140 à 340 Litres



# **AUTOMATIC AUTOCLAVE**

***Horizontal / Free Stand Type***

*Automatic Cycle Control by Electric Process Timers & Solenoid Valves  
from Heating Up » Sterilizing » Drying » Exhaust » Cycle End*

**Exist in 3 models: 143 - 196 - 339 LITERS**





# AUTOCLAVE ELECTRONIQUE

## 140 à 340 Litres

**Description:**

Type	Horizontal		
Voltage	AC 230/380/415V		
Control	PLC control		
Power supply	7KW	9KW	12KW
Designed temperature	140°C		
Sterilizing temperature	105~136°C		
Sterilizing pressure	0.2~2.3 kg/cm <sup>2</sup>		
Chamber size	45X90cm	50X100cm	60X120cm
Chamber capacity	143 Liters	196 Liters	339 Liters
Program	Wrapped/unwrapped/liquid		
Auto dry	Yes		
Auto add water	Yes		
Safety device	²Yes		
Vacuum drying	Yes		
Door interlock	Yes		
Printer	Yes		
N.W.	325 KGS	350 KGS	480 KGS
G.W.	425 KGS	500 KGS	630 KGS
Packing (LxWxH)	149X99X206cm	164X107X210cm	186X117X210cm

**Features:**

Heating element	1pc
Silicone door seal	1pc
SS304 tray	1pc
viewer tube	2pcs

**Safety Approval:**

