




AUTOCLAVE 25 LITRES



TABLE TYPE AUTOCLAVE

**Ideal for use in Medical clinics, Health centers,
Hospital departments, Dental clinics, Acupuncture
clinics, Laboratories...**



CE 0434
ISO 9001
ISO 13485



AUTOCLAVE 25 LITRES

SPECIFICATIONS	
MODEL NO.	25 L
CHAMBER SIZE	Stainless Steel #304 260mmØ x 450 mm (D)
CHAMBER CAPACITY	APPROX. 25 LITER
OVERALL	W540 X H450 X D630 mm
POWER CONSUMPTION	110V / 1.4 KW / 12.8 A 220V / 1.8 KW / 8.3 A
USING PRESSURE and TEMPERATURE	121°C / 1.2 KG/CM ² 134°C / 2.1 KG/CM ²
DESIGN TEMPERATURE	140°C
SAFETY DEVICES	PRESSURE CONTROL SWITCH DOOR PRESSURE AUTO-LOCK DOOR SAFETY CATCHER OVERHEAT PROTECTION SAFETY VALVE ELECTRONIC CIRCUIT SAFETY SYSTEM LOW WATER INDICATION AND CUT-OFF EMERGENCY EXHAUST SWITCH
SAFETY APPROVAL	CE-MDD by DNV CE-EMC, LVD by TUV CHAMBER DESIGNED BY ASME / U.S.A. STANDARD
STANDARD ACCESSORIES	Stainless Steel #304 Sterilization Box. / 1pc. Stainless Steel #304 Tray / 1pc.
OPTIONAL ACCESSORIES	Stainless Steel #304 Tray frame x 1 + Tray x 3pcs + Tray pulling holder x 1pc Stainless Steel #304 Sterilization spring holder x 2pcs.
AUTO-DRY FUNCTION	YES