



BUCKY STAND YZ-LX



Vertical Bucky Stand, which Bucky grid moved up and down, is a photograph device,
Moved up and down portable, manual lock freely
Support frame noiseless, cassette may direct into, operate convenience
Fixed Bucky no power (or electromotor oscillating grid)
Vertical movements 500mm~1600mm
Size of film 203mm × 254mm (8" × 10")
~356mm × 432mm (14" × 17")

Fixed Bucky X-Ray grid indexes

Film:	15" × 18"
40L/cm Density:	40L/cm
Grid proportion:	8:1
Converging focus distance fo:	1800mm