








CENTRIFUGEUSE 2060

■ CN-2060

- LED speed and time display.
- 10 memory programs for preset parameters.
- Microprocessor control.
- Auto accelerating and decelerating functions.
- Auto cover lock with motor running.
- Auto power off for rotor imbalance.
- Auto buzzer when finish.
- Max. speed: **6,000 rpm** ; Max. RCF: 3980 x g
- Auto power off when open.
- Auto-balance device.
- Max. capacity: 300 ml.
- Auto-cooling air circulation.
- ABS outer covering; see-through window.
- Other rotor selections for various tube size.



Rotor		Capacity (ml x pcs)	Max. RCF (x g)	Tube size (φ x l) mm	
				Max.	Min.
RA-1508		15 x 8	3,940	16.5 x 122	15 x 90
RA-1512S		15 x 12	3,940	16.5 x 122	15 x 90
RA-5004		50 x 4	3,980	29.5 x 115	27.5 x 90
RA-5006		50 x 6	3,980	29.5 x 115	27.5 x 90
REA-1520		1.5 x 20	3,140	11 x 48	10 x 40