



# CENTRIFUGEUSE HAUTE VITESSE

## ■ HIGH SPEED CENTRIFUGE CN-10000

- LED speed and time display.
- 10 memory programs for preset parameters.
- Microprocessor control.
- Auto power off for rotor imbalance.
- Auto buzzer when finish.
- High torque motor.
- Overload safety protection.
- Hidden auto-balance device.
- Auto power off when open.
- Auto cover lock with motor running.
- Power pilot light and fuse safety device.
- Auto accelerating and decelerating functions.
- Auto-cooling air circulation.
- See-through window.
- Optional rotor selections for various tube sizes.
- Low noise.



Model	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Timer (min:sec)	Weight (kg)	Size (cm) W x D x H
CN-10000	1,000 ~ 10,000	16,000	400	99:59	17	36 x 42 x 29

AC 110 / 220V, 50 / 60Hz



## CENTRIFUGEUSE HAUTE VITESSE

Applicable Model	Rotor		Max. Speed (rpm)	Max. RCF (x g)	Capacity (ml x pcs)	Tube Size (φ x L) mm		Cat. No.
						Max.	Min.	
CN-10000	RS-50		4,500	3,290	50 x 4TC	29.5 x 120	26.5 x 90	2-3400-23
					50 x 4	35 x 117	32 x 90	2-3400-22
			3,800	2,345	15 x 16	16.7 x 115	13.5 x 86	2-3400-32
								
	RA-1508		10,000	10,500	15 x 8	17 x 122	16 x 90	2-5052-31
	RA-1512H		9,000	10,130	15 x 12	17 x 122	16 x 90	2-3000-31
	RA-5004		8,000	6,720	50 x 4	29.5 x 117	27.5 x 90	2-5052-12
	RA-5006		7,000	5,400	50 x 6	29.5 x 117	27.5 x 90	2-5052-22
	RA-5008		6,500	5,280	50 x 8	29.5 x 117	27.5 x 90	2-1000-12