



CENTRIFUGEUSE HAUTE VITESSE

■ HIGH SPEED CENTRIFUGE CN-2200

- LED speed and time display; auto timer setting.
- 10 memory programs for preset parameters.
- Microprocessor control.
- Auto accelerating and decelerating functions.
- Auto cover lock when running and auto power off when open.
- Stainless steel rotor insert sterilizable.
- Max. capacity: 300 ml ;
- Timer: 99 min 59 sec.
- Auto buzzer when finish.
- Auto power off for rotor imbalance.
- ABS hard shell and see-through window.
- Auto-cooling air circulation.
- Wide rotors and accessories selections.



Model	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Timer (min:sec)	Weight (kg)	Size (cm) W x D x H
CN-2200	1,000 ~ 10,000	10,500	300	99:59	10	29 x 32 x 25

AC 110 / 220V, 50 / 60Hz

Model (Cat. No.)	Rotor		Max. speed (rpm)	Max. RCF (x g)	Capacity (ml x pcs)	Tube size (φ x L) mm		Cat. No.
						Max.	Min.	
CN-2200 (2-5053-03)	RA-5006		7,000	5,400	50 x 6	29.5 x 117	27.5 x 90	2-5052-22
	RA-5004		8,000	6,720	50 x 4	29.5 x 117	27.5 x 90	2-5052-12
	RA-1512S		9,000	8,900	15 x 12	17 x 122	16 x 90	2-5052-41
	RA-1508		10,000	10,500	15 x 8	17 x 122	16 x 90	2-5052-31



CENTRIFUGEUSE HAUTE VITESSE