



CENTRIFUGEUSE HEMATOCRITE

■ HEMATOCRIT CENTRIGUGE H-240

- Manual-braking system provides rapid deceleration.
- Regressive 15 min. timer with ring.
- Accessory: Hematocrit reader (order separately)
- Power pilot light and fuse safety protection.
- Auto-cooling air circulation.
- See-through window.
- Low noise.



Model	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Timer (min:sec)	Weight (kg)	Size (cm) W x D x H
H-240	12,000	1,880	15,000	0 - 15	8	29 x 32 x 25

AC 110 / 220V, 50 / 60Hz