



CENTRIFUGEUSE MICRO

■ MICRO MCD 2000

LED speed and time display.

Microprocessor control.

10 memory programs for preset parameters.

Auto accelerating and decelerating functions.

The transient centrifugal functions. Auto cover lock with motor running.

Auto power off for rotor imbalance.

Auto buzzer when finish.

ABS hard shell and see-through window.

Auto-cooling air circulation.

Low noise



Model	Speed (rpm)	Rotor	Max. RCF (x g)	Max. Capacity (ml x pcs)	Timer (min.sec)	Weight (Kg)	W x D x H Size (cm)
MCD-2000	1000-14,000	REA-1518	16,000	1.5/2.0 x 18	99:59	8	29 x 32 x 24.5