



CENTRIFUGEUSE pour IMMUNOLOGIE

■ IMMUNOLOGY BLOOD TYPING CENTRIFUGE CN-2160

- LED speed and time display.
- 10 memory programs for preset parameters.
- Auto-count RCF function.
- Effective RCF design: achieve required speed then begin countdown.
- Microprocessor control.
- Auto accelerating and decelerating functions.
- Auto buzzer when finish.
- Hidden auto-balance device.
- Auto cover lock with motor running.
- Power pilot light and fuse safety device.
- ABS outer covering; see-through window.
- Suitable for urine and blood samples.
- Auto-cooling air circulation.



Model (Cat. No.)	CN-2160 (2-5052-03)		
Rotor (Cat. No.)	SERO ROTOR (2-5050-11)		
Speed (No.)	1	2	3
Max. RCF (x g)	500	1,000	1,000
Timer (sec.)	60	15	60
Capacity	6ml x 12pcs of glass tube (φ12 x L75 mm)		
Application	The identification of blood type, red blood cell & lymph cell cleaning, the handle of blood coagulation, and the separation of lymph cell. To observe the reaction of blood coagulation, cultivate and separate cell, match blood by intersection, blood separation, and anti-serum globulin test.		