



CENTRIFUGEUSE REFRIGEREES

LCD, temperature, speed and RCF setting also display.
 Temperature range : 0°C~40°C.
 10 memory programs for preset parameters.
 Microprocessor control.
 Auto power off for rotor imbalance.
 Auto buzzer when finish.
 High function and high effective control system.
 Overheat and overload safety protection.
 Stainless steel inside bowl.
 Breakdown indication and fuse safety device.
 Auto accelerating and decelerating functions.
 Auto-cooling air circulation.
 Low noise.

LARGE SEIZE CN-650D

■ CN-650R



Model	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Timer (min:sec)	Weight (Kg)	W×D×H Size (cm)
CN-650R	500~6,500	4,900	1000	99:59	112	52×58×90

■ CN-650RD












LARGE SEIZE CN-650RD



Model	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Timer (min:sec)	Weight (Kg)	W×D×H Size (cm)
CN-650RD	500~6,500	4,900	1000	99:59	75	80×56×38



CENTRIFUGEUSE REFRIGEREE

Applicable Model	Rotor		Max. Speed (rpm)	Max. RCF (x g)	Capacity (ml x pcs)	Tube Size (φ x L) mm		
						Max.	Min.	
CN-650R	RS-250		3,500	2,250	15 x 96	10.7 x 100	8 x 41	
					6 x 112	12 x 110	8.5 x 60	
					10 x 80	15.8 x 100	12 x 60	
					15 x 40TC	16.5 x 122	14 x 60	
					15 x 48	17.2 x 122	14.5 x 60	
					50 x 16TC	29.5 x 125	27.5 x 50	
					250 x 4	62 x 130	58 x 70	
CN-650RD								
	RS-1548		3,500	2,430	15 x 48	17.5 x 107	14 x 80	
	RS-2504		3,500	2,480	250 x 4	65.5 x 130	61 x 90	
CN-650R	RS-9602		2,500	1,060	96 micro x 4	129 x 86	126 x 83	
	RS-100		3,500	2,090	6 x 56	13 x 105	10 x 50	
			3,500	2,150	10 x 48	16 x 107	13 x 50	
			3,500	2,150	15 x 16TC	16.7 x 135	13.5 x 86	
			3,500	2,150	15 x 32	17.5 x 115	14 x 80	
			4,000	2,770	50 x 4	35 x 135	32 x 90	
			5,000	4,640	50 x 4TC	29.5 x 137	26.5 x 90	
			3,500	2,120	50 x 8	35 x 120	32 x 90	
			3,500	2,120	50 x 8TC	29.5 x 127	26.5 x 90	
			4,000	2,930	100 x 4	45.5 x 130	42 x 95	
	CN-650RD							
	RS-50		4,500	3,280	15 x 16	16.7 x 115	13.5 x 86	
			4,000	2,600	50 x 4	35 x 117	32 x 90	
			5,500	4,900	50 x 4TC	29.5 x 120	26.5 x 90	
5,000			4,050	100 x 4	45 x 112	42 x 95		
								
RA-5008		6,000	4,500	50 x 8	29.5 x 118	27.5 x 90		
RA-1512H		6,000	4,500	15 x 12	16.7 x 130	13.5 x 86		