



CENTRIFUGEUSE REFRIGEREE

MEDIUM SEIZE

CN-350RD

- LCD, temperature, speed and RCF setting also display.
- Temperature range : 0°C~40°C.
- 10 memory programs for preset parameters.
- Microprocessor control.
- Auto power off for rotor imbalance.
- Auto buzzer when finish.
- High function and high effective control system.
- Overheat and overload safety protection.
- Stainless steel inside bowl.
- Breakdown indication and fuse safety device.
- Auto accelerating and decelerating functions.
- Auto-cooling air circulation.
- Low noise.














Model	Speed (rpm)	Max. RCF (x g)	Max. Capacity (ml)	Timer (min-sec)	Weight (Kg)	W X D X H Size (cm)
CN-350RD	500~6,500	4,500	400	99:59	52	65 X 53 X 34

Optional rotors and adapters:

See next page



CENTRIFUGEUSE REFRIGEREE

Applicable Model	Rotor		Max. Speed (rpm)	Max. RCF (x g)	Capacity (ml x pcs)	Tube Size (φ x L) mm	
						Max.	Min.
CN-350RD	RS-50		3,500	1,990	100 × 4	45 × 112	42 × 95
			4,000	2,600	50 × 4	35 × 117	32 × 90
			4,000	2,600	50 × 4TC	29.5 × 120	26.5 × 90
			3,500	1,990	15 × 16	16.7 × 115	13.5 × 86
			3,500	1,990	15 × 12TC	17 × 120	15 × 85
			3,500	1,990	10 × 24	15.5 × 113	14 × 85
			3,500	1,700	6 × 40	12.7 × 88	11 × 63
			3,500	1,700	6 × 36	12 × 106	11 × 85
			3,500	1,330	15 × 48	11 × 56	10 × 36
							
	RS-1508		3,500	1,990	15 × 8	17 × 122	15 × 90
	RA-1508		6,000	3,940	15 × 8	16.7 × 122	15 × 90
	RA-1512H		6,000	4,500	15 × 12	16.7 × 122	15 × 90
	RA-1516		4,000	2,250	15 × 16	16.7 × 115	15 × 90
	RA-1024		3,500	1,840	10 × 24	16.7 × 107	15 × 90
	RA-5004		6,000	3,980	50 × 4	29.5 × 130	27.5 × 94
	RA-5006		6,000	3,980	50 × 6	29.5 × 130	27.5 × 90
	RA-5008		6,000	4,500	50 × 8	29.5 × 117	27.5 × 90
	REA-1520		6,000	3,140	15 × 20	11 × 56	10 × 36