



CENTRIFUGEUSE REFRIGEREE

MICRO MCD-250RD

- LCD, temperature, speed and RCF setting also display.
- Temperature range : 0°C~40°C.
- 10 memory programs for preset parameters.
- Microprocessor control.
- Auto power off for rotor imbalance.
- Auto buzzer when finish.
- High function and high effective control system.
- Overheat and overload safety protection.
- Stainless steel inside bowl.
- Breakdown indication and fuse safety device.
- Auto accelerating and decelerating functions.
- Auto-cooling air circulation.
- Low noise.



| Model | Speed (rpm) | Rotor | Max. RCF (x g) | Max. Capacity (ml x pcs) | Timer (min:sec) | Weight (Kg) | W x D x H Size (cm) |
|-----------|-------------|----------|----------------|--------------------------|-----------------|-------------|---------------------|
| MCD-250RD | 1000~14,000 | REA-1518 | 16,000 | 1.5/2.0 x 18 | 99:59 | 35 | 55 x 42 x 28 |