



COAGULATION ANALYSER EXY

Photoelectricity-magnetolectric method, avoid interfering

Automatic electric sample-added (optional for EXY-20)

Open reagent system, it's more easy to find right reagent

Exist in 2 version: EXY-20 (2 Channels, EXY-40 (4 Channels)

Specifications

Testing Channel: 4 Channels (EXY-20: 2 Channels)

Analysis Parameters: PT, APTT, TT, FIB, Heparin, LMWH, PC, PS, D-Dimer, Factor II, Factor V, Factor VII, Factor VIII, Factor IX, Factor X, Factor XI, Factor XII

Sample Position: 12 (EXY-20: 6)

Reagent Position: 5

Printer: Low noise, long life, high speed thermal

Display: Large LCD

Data Communication: RS-232 serial cable

Data Storage Capacity: 99*4 test results, and power off protection

Power supply: 220VAC, 50-60Hz, 60W

Dimensions: 360 * 300 * 110mm

