



.DEFIBRILLATOR WITH SCOPE

DeltaLife 550 combines in a single unit the performances of an adult/pediatric synchronisable defibrillator, a monitor, ECG recorder and, optionally, a pacemaker and pulse oximetry.

DeltaLife 550 is a light and very reliable defibrillator, capable of 360 J maximum shock energy. The unit is ready to defibrillate in less than 10 seconds and can provide more than 50 shocks at maximum energy. A special version for internal defibrillation can be provided.



DeltaLife 550 is a simple and complete monitor. The high resolution LCD display shows clear ECG curves, heart rate value and all alarm status.

DeltaLife 550 is a digital ECG for 7 or 12 leads recordings. The unit provides both manual and automatic thermal printing of high resolution fully annotated ECG channels.

Technical specifications:

General

Dimensions	465 x 425 x 220 mm.
Weight	9 Kg
Power Supply	AC: 125/220 V. 50 or 60 Hz.
Battery	12 V Ni-Cd

Monitor - Display LCD

Dimensions	320 x 240 dot (115 x 86 mm.)
Visualization	4.5 sec. or 9 sec.
Sensitivity	5, 10, 20 mm/mV.
Visual and acoustic alarms	present
Autonomy	>3 hours

Recorder

Type	on 48 mm thermal paper
Speed	10, 25 and 50 mm/sec.
Operation mode	manual, automatic

Defibrillator

Energy	0 - 360 J
Charging time	10 seconds for 360 J
Available shocks	> 50 shocks at 360 J
Defibrillation paddles	adult and pediatric paddles (optional)