



LIT D'ACCOUCHEMENT

DELIVERY BED LM-01



Construction made of steel sections, coated with a powder varnish. Removable covers of seat segment and thigh segment. Adjustment of bed frame height is executed by foot-operated hydraulic pump (LM-01.0) or by electro-mechanical actuator (LM-01.1). Trendelenburg, reverse Trendelenburg positions and inclination angle of back rest are facilitated by gas springs with interlocking. Central wheels brake with one directional wheel.

TECHNICAL DATA LM-01:

Total bed length x width:	2180 x 980 mm
Height adjustment:	560 to 930 mm
Trendelenburg/reverse Trendelenburg position:	12°/5°
Back rest inclination:	70°