



DENTAL UNIT 638-1

Very good price – specifications relation

Standard Configurations:

Auto thermostatic rinsing water supply
 Auto cuspidor flush and cup filler
 External water & air purification
 Mute motor controlled chair (three memories)
 Rotatable ceramic cuspidor
 Saliva suction (high & low)
 Adjustable operation light
 Adjustable headrest
 Triple-purpose syringe (cold, hot water): 2pcs,
 Low voltage X-ray film viewer
 Multifunctional 4-hole pedal switch
 Pure water bottle

Technical Specifications:

Power supply voltage: AC220V/230V
 50Hz/60Hz (customisable),
 Motor voltage: DC24V,
 Air pressure: 0.5MPa - 0.8MPa,
 Water pressure: 0.2MPa - 0.4MPa,

Packing:

Wooden Cartoon: 139X90X109cm, 1.36cubic metres, G.W.:200kgs,



Options:

Set of instruments: High speed handpiece: 2pcs,
 Low speed handpiece: 1pc,



Wireless Endoscope System

1/4" Sony Super HAD CCD
 High resolution, High sensitivity, High performance
 Compact Size, Easy to use and focus
 2 or 3 keys function control, the new button can control: Previous image,
 next Image, and deleted image from the camera directly,
 The camera can with a small monitor on it, that dentist can see the picture from
 Display amazing clear images in two different viewing modes: Live and Freeze
 Effective Pixels (HxV): 811x 508



Scaler ultrasonic

Handpiece is cast from titanium;
 Excellent mechanic function;
 With excellent piez-ceramic;
 Handpiece has high & long lasting power & generates little heat;
 Cable is silica gel tube, soft & durable;
 Output power:3W~20W; Frequency : 30kHz±3kHz;



Curing light

Working voltage : AC100~240V 50/60H
 Output: DC-5V
 Wavelength:450-480mm
 Light power more than 1000mW/cm2



Dental compressor

30 L