



# EMERGENCY VENTILATOR

## Emergency portable ventilator P-10

### Applications:

Used for first aid, ambulance and patients transport within hospital with carry-bag and small cylinder or used for first aid and ICU with mobile table and humidifier. It can be used for pediatric and adult patients.



Breath Rate	10 15 20 25 30 35bpm
I:E	2:1 1:1 1:1.5 1:2 1:2.5 1:3
Tidal Volume	300 ~ 1000 ml
W×D×H	380×120×240mm
Weight	3kg
AC Voltage	110V~240V
DC Voltage	12V
Inlet Battery	12V , 2Ah