




# ENTERAL FEEDING PUMP CY-300

- Correcting the flow rate automatically.
- Regulating the flow rate conveniently.
- Presetting dose delivered.
- Indicating the total volume.
- Bionic designed.
- Fast infusion controlled.
- Audio and video alarms.
- With internal battery.
- Lock the flowrate when running for preventing accident.
- Indicating Power type (AC or DC).

Flow rate:	1 ml/h-300 ml/h (1 ml/h increments)
Preset dose volume	0-999 ml (5 ml increments)
Indicate total volume	0-999 ml
Flow error:	10%
Internal battery	9.6V±2V
Power	AC 220V ± 22V 50 Hz±1 Hz
Power input	13.5VA
Net Weight	1.6kg
Operation Temperature	5-40°C
Relative Humidity:	70%



**Alarms:** Flow error (Empty, OCCLUSION), Low battery, AC power error, Preset dose delivered finished, Reminder alarm.



*Designed Closely To Clinic,  
High Quality, High Security,  
The Best Choice For Feeding*



CY-B



CY-A

Specified designed feeding set, CY-B spike set and CY-A feeding bag