



# MAMOGRAPH 3000



## Specifications

High voltage generator	High frequency constant voltage: 40kHz
	Voltage rang: 22-35kVp, 1kV step
	Max. Current: 100mA
	Max. Output power: 5Kw
	Exposure condition control: Manual/Auto
X-ray tube	Dual focus: 0.1mm/0.3mm
	Anode heat capacity: 150KHU
	Anode type: Rotating Anode
C-arm	Rotating angle: +/- 145°
	SID: 60cm Fixed
	Anti-scatter grid: Carbon based fine plate
	Cassette size: 18x24cm
Dimensions	Gantry: 650mmx950mmx1860mm
	Console: 456mmx327mmx880mm
	C-arm vertical movement range: 630mm
	Overall weight: 250kg
Compressor	Manual/Auto: Auto (motor)
	Compression force/method: Max.18Kg/Auto, three speeds, flexible, fast pressure release
	Compression force display: LED
	Compressor paddle movement range: 268mm
Power	220VAC,30A



Le fabricant et nous même, nous réservons le droit de modifier le produit en vue d'une amélioration