



# MEDICAL REFRIGERATOR

## Features:

Used in hospitals, drugstore, epidemic prevention, medical research etc storage of biological and pharmaceutical products  
 YC-182B, YC-308B, YC-500B

Horizontal type adopts top open single door or double door (for capacity over 400 liters), easy open

Interior cabinet made of corrosive resistance pattern aluminum sheet or stainless steel provides smooth and clean plate surface, good effect on heat exchange and not easy corrosion

Extra thick insulation effectively reduce loss of cold air, energy saving, long service life

Control and display system for inside temperature on the front of the casing, provide adjustment and check inside temperature and operating status without opening the door, beautiful, practical and graceful

Microprocessor temperature control, digital temperature display can adjust and view inside temperature without opening the door, beautiful, practical and dignified.



Model	YC-182B	YC-308B	YC-500B
Type	Chest		
Capacity	182L	308L	500L
Inside temperature	0°C~8°C		
Power supply	~220V/50Hz		
Power consumption	0.7 kW·h/24h	0.9 kW·h/24h	1.3 kW·h/24h
Current input	0.7A	0.8A	1.1A
Prevention electricity degree	B type equipment		
To harmfully enters the fluid protection degree	Common equipment		
Working way	Interim work		
Climate class	N		
Fuse standard	6.3A		
Refrigerant/Charging	R134a/70g	R134a/120g	R134a/130g
External dimension WxDxH (mm)	1064x620x925	1125x710x925	1705x710x925
Packing dimension WxDxH (mm)	1128x636x990	1176x726x1049	1756x726x1049
Net/gross weight	60kg/65kg	70kg/75kg	80kg/85kg