



SURGERY MICROSCOPE LZ12

Be applicable to the **microsurgery surgical operations**, such as facial features section, **dentistry, orthopedics, microsurgery, ophthalmology and gynecology...etc.**

FEATURE:

1.PRIMARY MICROSCOPE

Model: double binocular

Manification 5X-25X

Magnification: motorized continuity

Diameter of viewing field 55-11mm

Large objective lens focus length
200,250,300mm

Illumination: a cold coaxial light

2.SECONDARY MICROSCOPE

Magnification 4X(viewing field diameter)

Swing angle:180°

3.STAND

Model: spring balanced

Up down range:500 mm

X-Y coupling shift field:40 mm

Swing angle of cross arm:240°

4.OPTIONAL ACCESSORIES:

Single binocular teaching microscope, CCD, camera, video

