



SURGICAL MICROSCOPE G3

Designed for micro surgery in otology, laryngology, neurosurgery, gynaecology, thoracicoperation and zootomy.

GX-3 is a fine surgical microscope with full performance.

It has comfortable eyepieces with large View-field and high-eye-point and continuous **electric zoom**.

It also has the **video-recorder** and permitting four persons to observe at the same time.

SPECIFICATIONS:

	GE-3
Magnification	6~30 X 4.2~20.9 X
Eye piece magnification	12.5 X
Field of view , diameter (mm)	7~36 10~51
Working distance (mm)	200, 300
Adjustment (mm)	350
Extension in length of arm (mm)	700
Lamp	24V 150W
Cold light	YES

STANDARD SET

	GE-3
Floor stand	X
Binocular tube	X
Lecture scope	X
Cold light supply	X
Foot panel	X
Camera	X
Flash unit	X

