



# SURGICAL MICROSCOPE SX

Cold light system

Designed for micro surgery in otology, laryngology, neurosurgery, gynaecology, thoracicooperation and zootomy.

Permitting TWO persons to observe at the same time.

<b>SXE-1</b>	<b>SXP-1CS</b>	<b>SXP-1CT</b>
		

### SPECIFICATIONS:

	<b>SXE-1</b>	<b>SXP-1CS</b>	<b>SXP-1CT</b>
Magnification	4~7.5 X 8~14.5 X	6 X, 10 X	6 X, 10 X
Eye piece magnification	10 X, 20 X	10 X, 16 X	10 X, 16 X
Field of view , diameter (mm)	21~38 14~26	16 25	16 25
Working distance (mm)	200, 250, 300, 375	200	200
Adjustment (mm)	310	600	260
Extension in length of arm (mm)	800	500	200
Rotation angle of binocular tubes	90° 30°	360°	360°
Lamp	24V 150W	8V 50W	24V 150W
Cold light	YES		
Cold light with 2 illuminator		YES	YES
Measuring eyepieces		YES	YES

### STANDARD SET

	<b>SXE-1</b>	<b>SXP-1CS</b>	<b>SXP-1CT</b>
Floor stand	X	X	X
Binocular tube	X	X	X
Lecture scope	X		
Cold light supply	X		
Measuring eyepiece		X	X