



BASIC X-RAY



P-300 40 kW
P-500 50 kW
High Frequency



The basic configuration comes with a floor-mounted tube stand		
The very slim profile of the floor rail reduces interference with personnel or stretchers moving through the room to a minimum		
Furthermore, the entire system does not require any wall or ceiling attachments, making for trouble free and easy installation.		
This series systems are equipped with HF generators with power level from 30 to 50kW, depending on your requirements in terms of application and workload.		
All the units have integrated APR (Anatomic Programs) to reduce parameter setup times and maximize your workflow.		
Main Unit	Three-phase	380 (±10%) V
	Capacity	≥40kW / 50 kW
	Frequency	50Hz ± 1Hz
	Resistance	≤0.2Ω
	Frequency	30kHz ± 1Hz

Le fabricant et nous même, nous réservons le droit de modifier le produit en vue d'une amélioration



BASIC X-RAY

	Max. output	320 mA / 100 kV / 200 mAs 500 mA / 100 kV / 250 mAs	
	Voltage	40~125 kV (45steps)	
	mAs	1.0~250 (47 steps)	
Table	◆◆	T17K	T17F
	Tabletop Height	750mm	700mm
	Tabletop longitudinal movement	550mm	550mm
	Tabletop transverse movement	230mm	230mm
	Wall stand longitudinal movement	1570mm	1570mm
	Column rotation	Nil	360°
	Tube rotating around Y axis	-20°~+120°	-20°~+120°
	Bucky device movement	810mm	550mm
	X-ray tube Vertical movement range	≥1150mm	≥1300mm

Le fabricant et nous même, nous réservons le droit de modifier le produit en vue d'une amélioration