



RESPIRATEUR PORTABLE

Portable ventilator P-100D

Applications:

Used for wake up room (PACU). (salle de réveil)

Performance Specification

Modes	CMV(IPPV), A/C(SIPPV) , SIMV ,	
SPONT , PEEP		
Tidal Volume	50~2000ml	
Breath Rate	1~100 bpm	
PEEP	0~2KPa	
O2%	45~100%	
Sigh	1~10 per 100 breath	
Peep	10-60 hpa	
Trigger Level	-20~0hpa	Over Pressure Relief
	<6kpa	



Alarms and Indicators

High and Low Pressure Alarm, lack of gas, Power alarm, Parameter error
Alarm Silence 60 sec

Monitors

Large LCD screen :

Data: Tidal Volume, O2%, Sigh, Respiratory Rate, Flow, Pressure (Synchronized Trigger, Peak Airway, CPAP, PEEP)

Curve: Airway Pressure Time Waveform Display

Electrical Specifications

AC Voltage: 110V~240V DC Voltage: 12V Inlet Battery:12V , 2Ah (Can be used for 8 >H)
Frequency: 50Hz/60Hz Power: 2A

Physical Specifications

Dimensions: Width: 300mm Depth: 160mm Height: 230mm Weight 3Kg