



ROLLATORS

Model 3 wheels	JM-01001A	JM-01001B
Overall width (mm)	660	660
Folded width (mm)	250	250
Height adjustment(mm)	800-920	800-920
Overall length (mm)	580	580
Weight limit for handle(kg)	120	120
Weight (kg)	6.7	6.7
Wheel	8"*2"TPR tyre	8"*2"TPR tyre
		

Model 4 wheels	JM-02001	JM-02002	JM-02003
Wide of front wheels (mm)	495	495	520
Height to handle(mm)	865-975	865-975	865-975
Seat height(mm)	570	570	570
Overall width(mm)	585	585	590
Overall length(mm)	590	590	590
Width between handles(mm)	470	470	460
Weight limit for handle(kg)	110	110	110
Weight(kg)	6	6	6
Wheel	6"*1" solid tyre	6"*1.2" solid tyre	6"*1.2" solid tyre
			



ROLLATORS

Model 6 wheels	JM-02004	JM-02005	JM-02006
Wide of front wheels (mm)	520	495	495
Height to handle(mm)	845-960	865-975	865-975
Seat height(mm)	570	570	570
Overall width(mm)	590	585	585
Overall length(mm)	590	590	590
Width between handles(mm)	460	470	470
Weight limit for handle(kg)	110	110	110
Weight(kg)	6	6	6
Wheel	5"*1" solid tyre	5"*1" solid tyre	5"*1" solid tyre
			

Model 6 wheels	JM-02011	JM-02012
Wide of front wheels (mm)	380	630
Height to handle(mm)	760-960	890-1000
Seat height(mm)	610	535
Overall width(mm)	610	460
Overall length(mm)	6800	710
Width between handles(mm)	435	670
Weight limit for handle(kg)	130	250
Weight(kg)	8	12
Wheel	8"*2"TPR tyre	7.5"solid tyre
		

Other models on request