



STERILISATEUR

STERILISATEUR – ETUVE – POUPINEL 60 litres

Digital display setting and actual temperature.
Microprocessor Control.
Double glass window.
Automatic alarm after operation.
Digital display running time.
Forced air circulation
PID self tuning, over temperature alarm
Stainless steel interior, with semicircle corners chamber.

Working size(mm) : 450x350x450
Power : 1000W
Temperature range : 50~300°C
Temperature resolution : 0.1°C
Temperature fluctuation : $\pm 1^\circ\text{C}$
Timer : 0-999min
Voltage : 110V 60Hz / 220V 50Hz
Dimension(mm) : 620x710x620
Packing info:
Size : 72cm*75cm*75cm
Gross Weight : 57kg
Net Weight : 48kg

