



SYRINGE PUMP A5803

Features:

- advanced technology and high safety standard
- wide range of flow rates
- automatic recognition of syringe size
- enhanced alarm functions
- internal rechargeable battery
- bright display
- simple and safe usage
- compact, light -weight and cost-effective



Ambient Temperature	5°C - 40°C
Power Supply	AC 220V 50Hz or DC 12V (AC110V available if required)
Power Consumption	Less than 25VA
Applicable syringes	20ml,50/60ml or 100ml disposable syringes; Compatible with most brands syringes
Infusion Limitation Options	20ml syringe 0.1ml/h-150ml/h; 50/60ml syringe 0.1ml/h-300ml/h; 100ml syringe 0.1ml/h-300ml/h; in 0.1 ml/h increments
Syringe Accuracy	Mechanical accuracy including syringe ± 5%
Alarms	Audible alarms with visual indicators of "NEAR EMPTY", "COMPLETE", "OCCLUSION", and "LOW BATT"
Included Accessories	a) 3-core standard power cable b) IV pole clamp c) DC 12v 2-pin plug d) Disposable Tubings
Battery Operation	4 hours continuously operation at 50 ml/hr with built-in rechargeable Ni-MH battery
Dimensions	L × B × H: 335×142×120mm
Weight	2.5 kgs / 5.5 lbs