



TABLE VETERINAIRE

- ◁ **Table Top** – Stainless steel made, and one drain hole on the right hand side corner of the table. 12 hooks under the side of table to tie up the pets during the surgery.
- ◁ **Base** – Rectangular, made of cast iron and cover with one piece of stainless steel shield. Base supported by four conductive rubber casters for mobility, and also with a locking screw to immobilize the table.
- ◁ **Table Top Control** – Controls of trendelenburg and reverse trendelenburg, lateral tilt.
- ◁ **Construction** – All exposed parts are either stainless steel, chromium plated or finished baked enamel for easy cleaning.



SPECIFICATIONS		
	Table Top Dimension	L123 X W60cm
	Height Adjustment	74~104cm
	Trendelenburg	28 ♦
	Reverse Trendelenburg	28 ♦
	Lateral Tilt	25 ♦
	Table Rotation	360 ♦
	Max Raise Weight	250 KGS
	N.W.	150 KGS
	G.W.	200 KGS
	Packing Size	L140 X W80 X H90 cm