



TREATMENT CHAIR 730



TREATMENT CHAIR

Rigid construction, chair is made of aluminium alloy, casting iron, outer case is made of fiberglass in match with aluminium alloy and coated with baked enamel.

Seat, back section and headrest made of soft rubber cushions.

Easy and simple operation, just use foot pedals for up and down function and manual control for back section reclining.



ISO 9001
ISO 13485


RTC-730

Designed For Better Life



TREATMENT CHAIR 730

SPECIFICATIONS

Chair Dimensions	L67 X W58 X H51~78 cm	Back Section Reclining	45 
Power Supply	AC 110/220V	N.W.	100 KGS
Electric Current	10/5A	G.W.	150 KGS
Chair Max. Elevation	27cm	Packing size	L87 X W77 X H135 cm