



VAPORIZER

For	Isoflurane	1208
	Enflurane	1208
	Halotane	1209

Capacity:Volume at Max mark:
Flow Range:
Temperature Range:
Concentration Range:
Maximum Pressure Load
Maximum Gradient:
Construction:

200ml nominal Volume at Min mark: 65ml
Operating flow range: 0.5~10L/min
Operating temperature range: 13~35C
0.5~5 vol%,relative error: +,- 20%
20.0kPa(Excess the Atmospheric Pressure)
45 angle
metal

