



VOLUMETRIC INFUSION PUMP A5808

AJ 5808 Volumetric Infusion Pump



Flow Rate Range	0.1-1200 ml/h
Flow Rate Increments	0.1ml/hr-10ml/h at 0.1 ml/h increments 10ml/hr-200ml/h at 1 ml/h increments 200ml/hr-1200ml/h at 10 ml/h increments
Total Volume Infused Display	0.1-9999ml
KVO Flow Rate	0.1ml/h ; 1ml/h ; 2ml/h
Purge Rate	1200ml/h
Accuracy	± 2% on pump only, ± 5% with IV set
Compatible I.V Set	Can be calibrated to fit all major brands
BOLUS	Adjustable Bolus rate from 0.1 to 1200ml/h with adjustable volume from 0.1 to 25ml
Display	Bright 2.7 inch OLED screen with separate large LED infusion rate display
Alarms, Audible and Visual	Infusion Completed, Door Open, Air-in-line, Occlusion, Empty Bag, Press Start Key, Battery Low, Battery Depleted, Initialization Error, Drive Mechanism Error
Multi-Language	English, French, Spanish
Power Source	AC 90V ~ 240V 50/60 Hz or DC 12V
Infusion Mode Memory	10 Preset infusion modes memory
Battery Capacity	Approximate 3 hours continuously operation at 100 ml/hr with built-in rechargeable NiMH battery pack (6 hours charging)
Net Weight	3.0 kg / 6.6 lbs
Dimensions	D × W × H: 172 × 105 × 240 mm D × W × H: 6 3/4" X 4" X 6 1/2"

Dual infusion mode



Adjustable Occlusion



Multi-language



Adjustable Bolus

