



DEFIBRILLATOR

Biphasic Defibrillator with monitor

BIPHASIC DEFIBRILLATOR MODEL LP-P

270 J

CE
0535



TECHNICAL SPECIFICATIONS

Wave Form	Biphasic	Paper Size	50mm x 30m
Operating Mode	Synchronous / Asynchronous	Max. Discharging	100+
Energy Selection	2 - 270 Joule	ECG Monitoring Time	4 Hours
ECG Derivations	DI-DII-DIII-AvR-AvL-AvF-C1...C6	Working Conditions	0/40 °C ve %30-95Rh
ECG Filters	Mon1-Mon2-Strong-Surgery-Diagnostic	Storing Conditions	-20 / +55 °C
Heart Rate Alarm Limits	15 - 300	Dimensions	403x152x324 mm
Freeze	Yes	Weight	8,5 Kg
Monitor Type	LCD	Class	II B
Monitor Resolution	320x240	Main Supply Voltage	AC230V/1A 50-60 Hz
Monitor Size	5.5 inch	Direct Working From Mains	Optional
Monitor Speed	25mm/sec	Pulse Oximeter	Optional
Printer Type	Thermal Printer	Standard Accessories	ECG Cable and Paper ECG Electrode - Power Cable
Printer Speed	25-50mm/sec		

Specifications & design subject to change without prior notice